

BERRY®

Love Life Always

Pulse Oximeter

BM1000B



Oxygen Level



Pulse Rate



Perfusion Index



Record in APP



Bluetooth



Shanghai Berry Electronic Tech Co., Ltd.

www.shberrymed.com

Tel: +86 21-58531958

E-mail: sales@berry-med.com

Unit 104, 1st Floor, 7th Building, NO. 1188 Lianhang Road, Minhang District, 201112 Shanghai, China.

BERRY®

Love Life Always

Pulse Oximeter


BM1000B



Technical Parameters

- SpO2 and pulse rate display.
- Perfusion index, pulse rate waveform and battery indicator.
- Real-time data transmission to PC and phone, real time storage in APP.

CE0123 NMPA ISO13485 FDA(510K)



BerrySmartHealth

Display:	TFT LCD	Working Voltage:	DC 2.2 V - DC 3.4V
Operation System:	IOS Android	Bluetooth:	BLE
Battery Type:	Two 1.5V AAA alkaline batteries	APP:	APP Store/Google Play: "BerrySmartHealth"
SpO2:	Detection scope: 35% - 100% Accuracy: ±2% 80% - 100% ±3% 70% - 79%	Pulse Rate:	Range: 25 - 250BPM Accuracy: ±2bpm Resolution: 1bpm
Size:	65(H)*32(W)*30(D)MM	Weight:	50g (include 2AAA alkaline batteries)
Operation environment:	Temperature: 5°C - 40°C	Relative humidity: 15% - 80%	Air pressure: 860hpa - 1060hpa

Shanghai Berry Electronic Tech Co., Ltd.

www.shberrymed.com Tel: +86 21-58531958 E-mail: sales@berry-med.com

Unit 104, 1st Floor, 7th Building, NO. 1188 Lianhang Road, Minhang District, 201112 Shanghai, China.